

OMRÄKNINGSTABELLER mmol/L, mg/dL

Omräkning mg/dL till mmol/L

(mg/dL x 0,0555 = mmol/L)

| mg/dL | mmol/L |
|-------|--------|
| 20 | ~ 1,1 |
| 40 | ~ 2,2 |
| 50 | ~ 2,8 |
| 60 | ~ 3,3 |
| 70 | ~ 3,9 |
| 80 | ~ 4,4 |
| 90 | ~ 5,0 |
| 100 | ~ 5,6 |
| 120 | ~ 6,7 |
| 140 | ~ 7,8 |
| 160 | ~ 8,9 |
| 180 | ~ 10,0 |
| 200 | ~ 11,1 |
| 220 | ~ 12,2 |
| 240 | ~ 13,3 |
| 260 | ~ 14,4 |
| 280 | ~ 15,5 |
| 300 | ~ 16,7 |
| 350 | ~ 19,4 |
| 400 | ~ 22,2 |
| 450 | ~ 25,0 |
| 500 | ~ 27,8 |

Omräkning mmol/L till mg/dL

(mmol/L x 18,016 = mg/dL)

| mmol/L | mg/dL |
|--------|-------|
| 2 | ~ 36 |
| 3 | ~ 54 |
| 4 | ~ 72 |
| 5 | ~ 90 |
| 6 | ~ 108 |
| 7 | ~ 126 |
| 8 | ~ 144 |
| 9 | ~ 162 |
| 10 | ~ 180 |
| 11 | ~ 198 |
| 12 | ~ 216 |
| 13 | ~ 234 |
| 14 | ~ 252 |
| 15 | ~ 270 |
| 16 | ~ 288 |
| 17 | ~ 306 |
| 18 | ~ 324 |
| 19 | ~ 342 |
| 20 | ~ 360 |
| 25 | ~ 450 |



Gilla GlucoMen på Facebook
www.areo.nu
shop.menarinidiagnostics.se

