

**A.MENARINI**  
diagnostics

www.menarindiagnostics.se

**SNABBA  
OCH PÅLITLIGA  
RESULTAT!**



**Produktkod (25 testkit):  
51674 COVID-19 Snabbtest-kassett**

**Endast 150 kr/test  
3750 kr/för 25 test (1 fp)**

**OBS endast till vården**

**Beställ: info@menarindiagnostics.se**

## COVID-19 IgG/IgM SNABBTEST


COVID-19 IgG/IgM snabbtest – ett kombinerat IgM/IgG antikroppstest för Covid-19 som används för att snabbt och kvalitativt upptäcka förekomst av IgG- och IgM-antikroppar av det nya coronaviruset. Validerat i ett flertal studier, bl a från Uppsala och England<sup>\*(1)</sup>

### Funktioner och Fördelar

- Snabba och pålitliga resultat (på endast 2–10 minuter)
- Enkelt, tidsbesparande förfarande - ingen specialutbildning behövs
- Test kan utföras med helblod (venös och kapillär), serum eller plasma
- Liten provstorlek: endast 5 µL av serum/plasma eller 10 µL helblod
- Alla nödvändiga reagenser ingår och ingen utrustning behövs
- Hög sensitivitet (90,11%) och specificitet (98,46%)
- Testet kan lagras vid rumstemperatur

**COVID-19 IgG/IgM Rapid Test Cassette utvärderades kliniskt av Renmin Hospital - Wuhan University, Kina.<sup>\*(2)</sup>**

## Testförfarande och tolkning



<sup>\*(1)</sup> Hoffman et al: Evaluation of a COVID-19 IgM and IgG rapid test; an efficient tool for assessment of past exposure to SARS-CoV-2  
Pallet et al: Point-of-care serological assays for SARS-CoV-2 in a UK hospital population: potential for enhanced case finding

<sup>\*(2)</sup> Clinical Evaluation Report of COVID-19 IgG/IgM Rapid Test Cassette (Whole Blood/Serum/Plasma)